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The biological compatibility of intelligent textiles

Von Dr. Corinna Petsch, Head of BioLabs, senetics healthcare group GmbH & Co. KG
By Dr Corinna Petsch, Head of BioLabs, senetics healthcare group GmbH & Co. KG

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SMART FLEXIBLES – COMFORTABLE HEALTHMONITORING OF THE FUTURE.

SmartFlexibles are at the latest since the development of the first flexible smartphone and are discussed and tested for numerous applications. In addition to the flexible smartphones, there are also textiles that can detect various parameters of the human body and enable continuous monitoring of the pulse rate, blood pressure and body temperature. This type of textiles is gaining in importance, especially in the gym sector and the health care.

CONTINUOUS SKIN CONTACT INVOLVES RISKS. However, since these materials are in direct contact with the skin for the most time, the question of biological compatibility is extremely important. In addition to the direct contact, transpiration also carries risks during the application. Due to its composition, human perspiration can attack the textiles and remove critical substances from the tissue. These can easily penetrate into the body through opened pores and the direct skin contact, which might cause undesirable side effects.

The effects of poor biocompatibility can be varied and depend on the field of application. For example, skin rashes and allergies can occur during skin contact. These are unpleasant for the user and lead to the product being no longer being used. This can have serious consequences for devices used for health monitoring.

WHAT IS BIOCOMPATIBILITY? Biocompatible materials and substances do not have a demonstrable negative effect on organisms and living beings in their environment. Biocompatible materials can also be divided into three categories: biotolerant, bioinert and bioactive (see Fig. 1).

Biotolerant products are best tolerated. They can remain in the body for years without being toxic. In contrast, interactions can be observed with bioinert substances. These are, however, so mild that no limits are exceeded at any time. Also, no rejection reactions can be detected. The term bioactive applies primarily to

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
Fig. 1: Definition of biocompatibility

Biocompatibility can be divided into three sub-areas. Biotolerant products can remain permanently in the body without toxic effects. In contrast, low interactions within the limits can be observed for bioinert substances. Bioactive materials interact with tissues and bones.

implants that interact with bone. These can be coatings which allow the bone to adhere to the implant. All other materials which have an influence on the organism are called biologically incompatible. In this case, clear foreign-body reactions can be observed in the organism. The effects of intolerance differ in their severity and depending on the field of application.

PROOF OF BIOLOGICAL COMPATIBILITY. Various in-vitro and in-vivo tests are available for the detection of biological compatibility. These are specified by the international ISO series 10993. It is valid for medical products, biomaterials and materials which have direct or indirect contact with the human body. The norm covers 20 sub-chapters and takes account of regulatory-organizational aspects as well as biological, chemical and physical tests. The type and number of tests dependent on the final use of the product. The following aspects should be taken into account when selecting the tests: the type of contact, the duration of contact and the risk potential (Fig. 2). For the type of contact, a distinction is made between superficial (skin) and internal (internal organs) contact and intrusion from outside into the body. The contact time must


Fig. 2: Aspects for selecting the necessary tests

Fundamental aspects to be considered during the selection are the type of contact, the duration of contact and the risk potential. The number and the effort of the tests differ and depend on the final use.

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also be considered accordingly. Thus, implantable products require more numerous and more severe tests due to the field of application. In addition, a detailed risk assessment is necessary in advance. This includes, among other things, a detailed literature review with attention to similar, already established materials and previous tests.

The following three tests are recommended for textiles which are in direct contact with the skin, but do not have contact with the body's interior: tests for in-vitro cytotoxicity (DIN EN ISO 10993-5) and tests for irritation and skin sensitization (DIN EN ISO 10993-10).

In principle, the question of the cytotoxicity of the material should first be clarified. This test is checked with established cell lines in a viability assay (see Fig. 3). The material is rinsed in medium in order to wash out soluble constituents, which could also dissolve during later use. This extract is then pipetted onto murine cells in various dilution stages and remains there for 24 hours. The evaluation is then carried out on the basis of the metabolic activity of the cells. If a high metabolic turnover is observed, the cells are vital and no cytotoxic effect can be detected. If there is a cytotoxic effect, the cells die and thus a reduced or no metabolic activity can be detected.


Fig.3: MTT cytotoxicity test with murine cells
In the case of positive control, all cells are dead, so that no color change and therefore metabolism can be detected. In contrast, the cells of the negative control are viable and metabolism can be shown. On the left and on the right are blind samples without cells; in the middle the extract of the test object was applied in various concentrations. The cytotoxicity test shown is negative, cell metabolism is detectable and therefore the extract shows no cytotoxicity potential.

If it is found that the material is not cytotoxic, further tests can be carried out with regard to skin irritation and sensitization. For this purpose, the effect of the extracts on 3D skin models is examined and the lymphocyte proliferation in murine lymph nodes is determined. If there are signs of negative interactions, the material selection should be reconsidered. The conduction of these tests is decisive for the success of the product, since this is the only way to ensure an irritant-free wearing of the textiles.

BIOCOMPATIBILITY ESSENTIAL FOR USER SAFETY. In summary, it can be stated that biocompatibility tests should be taken in advance to avoid incompatibilities. This saves costs and time, which require

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senetics healthcare group GmbH & Co. KG, D-91522 Ansbach,
www.senetics.de

possible material changes. If, during the course of the investigations, it becomes apparent that a raw material or material is not compatible, this can have an effect on the entire development process and, in the worst case, require a complete restart of the project!

The testing of biological compatibility should therefore be regarded as an essential component of product development and should not be taken lightly! Careful planning and testing during development can save valuable time and unnecessary costs.

senetics healthcare group GmbH & Co. KG, D-91522 Ansbach,
www.senetics.de